



# SURGICAL GOWNS

# SURGICAL GOWNS



## AERO CHROME\* BREATHABLE PERFORMANCE SURGICAL GOWN

AAMI Level 4

### Viral Protection Gown

- AAMI Level 4 Liquid Barrier Protection<sup>1</sup>
- Viral Protection ASTM F1671<sup>2</sup>
- EN 13795 High Performance<sup>1</sup>
- Clearly identifiable Silver color
- Soft, Light and Breathable SFSMS Fabric<sup>1</sup>
- Raglan Sleeves
- 2 scrim hand towels
- Tested with 52 Chemotherapy Drugs and MABs<sup>3</sup>

Code	Description	Packaging
<b>48395</b>	Aero Chrome* Breathable High Performance Surgical Gown with Towels, L	Items/Case = 30
<b>48396</b>	Aero Chrome* Breathable High Performance Surgical Gown with Towels, XL	Items/Case = 28
<b>48397</b>	Aero Chrome* Breathable High Performance Surgical Gown with Towels, L, X-Long	Items/Case = 30
<b>44689</b>	Aero Chrome* Breathable High Performance Surgical Gown with Towels, L	Items/Case = 30
<b>44690</b>	Aero Chrome* Breathable High Performance Surgical Gown with Towels, XL	Items/Case = 28
<b>44691</b>	Aero Chrome* Breathable High Performance Surgical Gown with Towels, XXL	Items/Case = 28



## AERO BLUE\*

AAMI Level 3

### AERO BLUE\* Surgical Gowns with Towels

- AAMI Level 3 Liquid Barrier Protection<sup>4</sup>
- Passes ASTM F1670 -Synthetic Blood Penetration test<sup>5</sup>
- Lightweight Breathable Film laminate Fabric<sup>6,7</sup>
- Raglan Sleeves
- Colour-Coded Neck Binding
- Book Fold

Code	Description	Packaging
<b>41928</b>	AERO BLUE* , Sterile, Large	Items/Case = 30
<b>41929</b>	AERO BLUE* , Sterile, X-Large	Items/Case = 28
<b>41930</b>	AERO BLUE* , Sterile, Large, X-Long	Items/Case = 30
<b>41931</b>	AERO BLUE* , Sterile, X-Large, X-Long	Items/Case = 30
<b>41733</b>	AERO BLUE* , Sterile, Large	Items/Case = 32
<b>41734</b>	AERO BLUE* , Sterile, X-Large	Items/Case = 30
<b>41724</b>	AERO BLUE* , Sterile, Large, X-Long	Items/Case = 30
<b>41726</b>	AERO BLUE* , Sterile, X Large, X-Long	Items/Case = 30

All codes in **BOLD** indicate a sterile product.

**NOTE:** Please see package inserts for directions for use for each of these products.

# SURGICAL GOWNS




## VALUESELECT\*

AAMI Level 2

### VALUESELECT\* Surgical Gown, AAMI 2, with Towels

- AAMI PB: 70 AAMI Level 2 fluid barrier protection<sup>8</sup>
- Sleeve and shoulders are ultrasonically bonded for fluid protection<sup>9</sup>
- Lightweight gown of SMMMS fabric<sup>10</sup>
- 360-degree AAMI 2 protection, sterile surgical gown<sup>8</sup>
- Alcohol-repellant and Anti-static<sup>11</sup>
- Low linting as per EN ISO 9073-10 Lint test<sup>12</sup>
- Low abrasion as per test for abrasion ASTM D 4966<sup>13</sup>
- Meets Class 1 standard for Flammability test CPSC 16 CFR Part 1610<sup>14</sup>
- Passes skin Sensitization and Skin irritation test: ISO 10993-10-2010<sup>15</sup>
- Passes cytotoxicity test ISO 10993-5-2009<sup>16</sup>
- Set-in Sleeves
- 2 Scrim-reinforced lint-free hand towels<sup>9</sup>

Code	Description	Packaging
<b>77010</b>	ValueSelect Surgical Gown, Sterile, L	Items/Case = 50
<b>77011</b>	ValueSelect Surgical Gown, Sterile, XL	Items/Case = 50



## STANDARD/EVOLUTION\* 4

### EVOLUTION\* 4 Set-In-Sleeve Non-Reinforced Surgical Gown with Towels

- EVOLUTION\* 4 Fabric
- Set-In-Sleeves
- Colour-Coded Neck Binding
- Adjustable Hook & Loop Neckline
- Book Fold

Code	Description	Packaging
<b>78940</b>	Standard Gown M	Items/Case = 60
<b>78945</b>	Standard Gown L	Items/Case = 60
<b>78950</b>	Standard Gown XL	Items/Case = 60

References: 1-16 Data on file

All codes in **BOLD** indicate a sterile product.

**NOTE:** Please see package inserts for directions for use for each of these products.



**O&M Halyard India Pvt Ltd**  
World Trade Center, Tower-1 Unit 308  
Kharadi, Pune 411 014  
Maharashtra, India

For feedback/suggestions/complaints, contact Customer care officer  
Tel no: +91 20-46845500/5555  
Email: [customercare\\_IND@hyh.com](mailto:customercare_IND@hyh.com)

COPY-04214

O&M Halyard makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

\*Registered Trademark or Trademark of O&M Halyard Inc. or its affiliates.  
©2021. All rights reserved

[www.halyardhealth.in](http://www.halyardhealth.in)  
 **O&M Halyard APAC**